


 quattro **DR^{HD}**



A whole new category of Digital Radiography...

Cuatro DR HD (High Definition) Digital Radiography is now available exclusively from Heska. Cuattro DR HD delivers detail, resolution, and contrast image quality that is visibly better, much better, than Standard Definition systems.

 **cuattro DR^{HD}**



HD Minimum Standards:

- Direct Deposition Cesium
- 100 micron pixel size
- HD pixel density. 18 Million (17x17)
- HD pixel density. 15 Million (14x17)
- Full Size Image Area (14x17 - 17x17)
- Resolution of at least 5 lp/mm
- 16 Bit A/D 16 Bit Grayscale Native
- 16 Bit Multi-Frequency Tuning



Cuattro DR HD Added Benefits:

- 24/7 Live Technical Support
- Automatic Exposure Detection
- Cassette Sized Form Factor
- Embedded - 64 bit OS
- Multi-Touch User Interface
- Cloud Based PACS Connectivity
- Cloud Archival
- Software Upgrades and Updates

A smarter direction in Digital Radiography. Better, much better.

www.heska.com.au • 1300 437 522 • imaging the possibilities

© 2020 Heska Corporation. All rights reserved. HESKA is a registered trademark and Cuattro DR HD is a trademark of Heska Corporation in the United States and in other countries. Specifications are subject to change without notice and may vary based on operating conditions. Pricing dependent on final configuration. AU20LT0103

