

Complete Test Menu

Element DC®, DCX™ and DC5X® Veterinary Chemistry Analysers

INDIVIDUAL SLIDES

BUN, CREA, ALP, ALT, PHOS, GGT, TBIL,
AST, GLU, TP, ALB, Ca, Mg, NH₃, TRIG,
LIP, AMY, CK, URIC, CHOL, TCO₂, LDH,
Electrolytes (Na⁺, K⁺, Cl⁻)

PANELS

PRE-SURGICAL PANEL (eWrap)

ALP, ALT, BUN, CREA, GLU, TP

NEW!

COMPREHENSIVE+ PANEL (eWrap)

ALP, ALT, BUN, CREA, GLU, TP, TBIL, ALB,
PHOS, Ca, Chol, GGT, LIP, AMY, Mg, TRIG,
AST, Electrolytes (Na⁺, K⁺, Cl⁻)

COMPREHENSIVE PANEL (eWrap)

ALP, ALT, BUN, CREA, GLU, TP, TBIL, ALB,
PHOS, Ca, CHOL, GGT

PLUS PANEL (eWrap)

LIP, AMY, Mg, TRIG, AST,
Electrolytes (Na⁺, K⁺, Cl⁻)

KIDNEY PANEL

ALB, BUN, CREA, Ca, PHOS, TP

LIVER PANEL

ALB, ALP, ALT, GGT, GLU, TBIL

EQUINE PANEL

ALB, AST, BUN, CREA, Ca, CK, GGT, GLU,
LDH, PHOS, TBIL, TP

CALCULATED VALUES

Globulin
Na⁺/K⁺ Ratio
Albumin/Globulin Ratio
Corrected Calcium
Estimated Osmolality
BUN/Creatinine Ratio

Element HT5® Veterinary Haematology Analyser

WHITE BLOOD CELL

Total White Blood Cell (WBC)
Neutrophil Number (NEU#)
Neutrophil Percent (NEU%)
Lymphocyte Number (LYM#)
Lymphocyte Percent (LYM%)
Monocyte Number (MON#)
Monocyte Percent (MON%)
Eosinophil Number (EOS#)
Eosinophil Percent (EOS%)
Basophil Number (BAS#)
Basophil Percent (BAS%)

RED BLOOD CELL

Total Red Blood Cell (RBC)
Haemoglobin (HGB)
Haematocrit (HCT)
Mean Cell Volume (MCV)
Mean Cell Haemoglobin (MCH)
Mean Cell Haemoglobin Conc. (MCHC)
Red Blood Cell Distribution (RDW-CV)

PLT

Total Platelet (PLT)
Mean Platelet Volume (MPV)

VISUAL RESULTS

WBC Distribution Scatterplots (3)
WBC Distribution Histogram
RBC Size Distribution Histogram
PLT Size Distribution Histogram

Element i+® Immunodiagnostic Analyser

MEASURED PARAMETERS

Total T4, TSH (Canine), Cortisol and
Progesterone

COMING SOON:

SDMA, Infectious Disease, cCRP, and more.

HeskaView Telecytology

SAMPLE TYPES

CBC Pathology, Fine Needle Aspirates,
Body Cavity Effusions, Joint Taps,
Bronchoalveolar Lavage, Transtracheal
Wash, Reproductive Cytology, and more

COMING SOON:

Element AIM™ Automated Artificial Intelligence Microscopy

URINE SEDIMENT ANALYSIS

RBCs, WBCs, Bacteria, Crystals, Casts,
and more

FECAL FLOTATION ANALYSIS

Roundworm, Hookworm, Tapeworm,
Whipworm, Giardia, and more



For further assistance, please call Heska's Technical Support Services

US 800 464 3752
www.heska.com

CA 866 382 6937
www.heskavet.ca

AU 1300 437 522
www.heska.com.au