

NEW



trūRapid Series

Infectious and Vector Borne Disease Testing

Introducing Heska® trūRapid™, the latest point-of-care infectious disease testing for all practices.

Focus on prevention and early detection with fast, accurate, and easy-to-use trūRapid tests. Bloodborne and faecal infectious disease diagnostics - now in the palm of your hand.



10 minutes to results

Room temperature storage

Excellent specificity and sensitivity

Long shelf life

www.heska.com.au



trūRapid Series

Rapid Detection In-house

...with accurate results in 10 minutes.

The Heska trūRapid line offers veterinary professionals an ease and breadth of testing, which includes as many as 8 common companion animal infectious diseases – all returning results in around 10 minutes. With the additional advantages of room temperature storage and long shelf-life, Heska trūRapid tests are versatile, reliable, and cost-effective.

Utilise multiplexed options like trūRapid FIV/FelV and trūRapid Parvo/Giardia/Coronavirus for more speed and flexibility!



| Test Name | Detection | | |
|---|--|---|-------------------|
| trūRapid Heartworm (Dirofilaria) | <i>Dirofilaria immitis</i> antigen | ■ | Canine |
| trūRapid Distemper Ab | Canine Distemper Virus antibody | ■ | |
| trūRapid FIV/FelV | Feline Leukemia Virus antigen / Feline Immunodeficiency Virus antibody | ■ | Feline |
| trūRapid FIV | Feline Immunodeficiency Virus antibody | ■ | |
| trūRapid FelV | Feline Leukemia Virus antibody | ■ | |
| trūRapid Feline Titer Test * | Feline Panleukopenia Virus antibody | ■ | Canine/ Feline |
| trūRapid Parvovirus Giardia Coronavirus | Canine Parvovirus / Coronavirus / <i>Giardia duodenalis</i> antigen | ■ | |
| trūRapid Parvovirus | Canine Parvovirus antigen | ■ | |
| trūRapid Giardia | <i>Giardia duodenalis</i> antigen | ■ | |
| trūRapid Parvovirus Giardia | Canine Parvovirus / <i>Giardia duodenalis</i> antigen | ■ | |

Sample Type: ■ Whole Blood, Serum, or Plasma ■ Faeces

* Coming soon.