

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

ME0126A

2026-03-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0126A-1	LOT	ME0126A-2	LOT	ME0126A-3
element HT5+ QC Mode	Ret# ×10 ³ /μL	76.4	± 64.4	233.3	± 99.0	189.5	± 99.0
	Ret%	1.29	± 1.12	4.89	± 2.00	8.31	± 3.50
	IRF %	12.3	± 15.0	12.7	± 13.0	9.3	± 25.0
	LFR %	87.7	± 15.0	87.3	± 15.0	90.7	± 25.0
	MFR %	10.7	± 12.0	11.4	± 10.0	6.4	± 22.0
	HFR %	1.6	± 10.0	1.3	± 10.0	2.9	± 10.0
	RHE pg	24.3	± 12.0	22.6	± 12.0	19.8	± 12.0
	IPF %	7.5	± 15.0	7.5	± 15.0	7.6	± 15.0

Before using, refer to the instruction sheet for mixing directions.