

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0326A

 2026-05-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0326A-1	LOT	ME0326A-2	LOT	ME0326A-3
element HT5+	Ret# x10 ³ /μL	74.6	± 64.4	233.8	± 99.0	181.4	± 99.0
QC Mode	Ret%	1.25	± 1.12	4.87	± 2.00	8.10	± 3.50
	IRF %	12.0	± 15.0	12.6	± 13.0	8.9	± 25.0
	LFR %	88.0	± 15.0	87.4	± 15.0	91.1	± 25.0
	MFR %	10.4	± 12.0	11.1	± 10.0	6.1	± 22.0
	HFR %	1.6	± 10.0	1.5	± 10.0	2.8	± 10.0
	RHE pg	24.2	± 12.0	22.4	± 12.0	20.6	± 12.0
	IPF %	7.5	± 15.0	7.6	± 15.0	7.6	± 15.0

Before using, refer to the instruction sheet for mixing directions.