

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0726A



2026-09-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0726A-1	LOT	ME0726A-2	LOT	ME0726A-3
element HT5+	Ret# $\times 10^3/\mu\text{L}$	70.2	\pm 64.4	216.7	\pm 99.0	181.7	\pm 99.0
QC Mode	Ret%	1.19	\pm 1.12	4.64	\pm 2.00	8.22	\pm 3.50
	IRF %	12.3	\pm 15.0	11.3	\pm 13.0	8.1	\pm 25.0
	LFR %	87.7	\pm 15.0	88.7	\pm 15.0	91.9	\pm 25.0
	MFR %	10.9	\pm 12.0	9.3	\pm 10.0	4.9	\pm 22.0
	HFR %	1.4	\pm 10.0	2.0	\pm 10.0	3.2	\pm 10.0
	RHE pg	24.8	\pm 12.0	23.7	\pm 12.0	21.5	\pm 12.0
	IPF %	7.5	\pm 15.0	7.7	\pm 15.0	7.7	\pm 15.0

Before using, refer to the instruction sheet for mixing directions.