

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

ME0925A

2025-11-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0925A-1	LOT	ME0925A-2	LOT	ME0925A-3
element HT5+ QC Mode	Ret# ×10 ³ /μL	78.8	± 64.4	242.6	± 99.0	180.6	± 99.0
	Ret%	1.34	± 1.12	5.24	± 2.00	7.99	± 3.50
	IRF %	11.9	± 15.0	12.6	± 13.0	10.4	± 25.0
	LFR %	88.1	± 15.0	87.4	± 15.0	89.6	± 25.0
	MFR %	10.3	± 12.0	10.8	± 10.0	8.2	± 22.0
	HFR %	1.6	± 10.0	1.8	± 10.0	2.2	± 10.0
	RHE pg	26.1	± 12.0	22.6	± 12.0	20.7	± 12.0
	IPF %	7.8	± 15.0	7.6	± 15.0	7.6	± 15.0

Before using, refer to the instruction sheet for mixing directions.